



Technical Data

Tracoe Vario XL

Size	ID mm	OD mm	OB mm	CD mm	BA °	REF 451 / REF 451-P	REF 461 / REF 461-P	REF 471 / REF 471-P
						TL (min-max) mm	TL (min-max) mm	TL (min-max) mm
07	7.0	9.7	110	23.0	100	65.0 - 98.0	64.5 - 97.5	60.0 - 93.0
08	8.0	11.2	128	28.0	100	76.0 - 116.0	73.5 - 113.5	73.0 - 113.0
09	9.0	12.3	134	30.0	100	73.5 - 118.5	78.0 - 123.0	73.5 - 118.5
10	10.0	13.7	146	32.0	100	85.0 - 130.0	84.5 - 129.5	80.0 - 125.0

ID: inside diameter (clear width) at bottom of cannula; **OD:** outside diameter at bottom of cannula; **TL min:** minimal length along centre line; **TL max:** maximal length along centre line; **OB:** length along outer bend; **BA:** bending angle; **CD:** cuff diameter